

Floor Standing • Low Speed • Refrigerated • Centriuge

FR6-B6



Scan to watch product video



6000 rpm	Maximum Speed
6880 xg	Maximum RCF
6x1200 ml	Maximum Capacity
6 pcs 500ml blood bags	
≤58dB (A)	Noise Level(dBA)
-20℃~40℃	Temp.range Accuracy of ±1℃

- * RTC (Real-Time Control) Function
- * Castor wheels for ease of mobility

Technical Specifications:

Speed accuracy	±20r/min	Daily-use program	30
Timer range	0~99h 59min adjustable	Acceleration/ Deceleration rates	12/12
Drive System	AC variable frequency motor drive	Temp accuracy	±1℃
Dimension	840×730×1220mm(L×W×H)	Display	LCD
Power supply	AC220V, 50Hz, 30A, 2/1.5kw	Net Weight	260kg

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Fixed Angle Rotor	6×500ml	6000r/min	6880xg	PP/PC tube
NO.2 	Swing Rotor	6×1200ml	4200r/min	5600xg	Oval cup
NO.3 	Swing Rotor	6×1000ml	4200r/min	5600xg	Round cup
NO.4 	Swing Rotor	4×1500ml	4200r/min	5600xg	Oval cup