

Floor Standing • Low Speed • Refrigerated • Centriuge

FR8-B4



Scan to watch product video



7000 rpm	Maximum Speed
8897 xg	Maximum RCF
4x1000 ml	Maximum Capacity Swing Bucket
4x500 ml	Fixed Angle Rotor
≤58dB (A)	Noise Level(dBA)
-20℃~40℃	Temp.range Accuracy of ±1℃

- * RTC (Real-Time Control) Function
- * Castor wheels for ease of mobility

Technical Specifications:

Speed accuracy	±20r/min	Acceleration/Deceleration rates	1--10
Time range	1-99min59s adjustable	Daily-use program	30
Dimension	690×590×835(L×W×H)	Net weight(without rotor)	190kg
Display	Touch LCD screen	Power supply	AC220V,50HZ,16A,1500W
Control&drive	High precision and low noise AC frequency conversion motor direct drive		
Microprocessor	Microprocessor control		

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Angle rotor	4×500ml	7000r/min	8897xg	PP/PC tube
	Adapter	4×100ml			PP/PC tube
	Adapter	4×3×50ml			PP/PC tube
	Adapter	4×8×15/10ml			PP/PC tube
NO.2 	Angle rotor	16×50ml	7000r/min	7121xg	PP/PC tube
	Adapter	16×15ml			PP/PC tube
	Adapter	16×10ml			PP/PC tube
	Adapter	16×3×1.5ml			PP/PC tube
NO.3 	Angle rotor	32×15ml	7000r/min	7300xg	PP/PC tube
NO.4 	Swing bucket	4×1000ml	4000r/min	3578xg	PP/PC tube
	Adapter	4×26×5/2ml			PP/PC tube
	Adapter	4×13×15/10ml			PP/PC tube
	Adapter	4×5×50ml			PP/PC tube
	Adapter	4×3×100ml			PP/PC tube
NO.5 	Microplate rotor	4×2×96well plate	4200r/min	2859xg	PP/PC tube
	Adapter	4×32×1.5ml			PP/PC tube