

Gel Card Centrifuge

MK24



Scan to watch product video



4000 rpm	Maximum Speed
2773 xg	Maximum RCF
4x12 cards / box	Capacity
≤58dB (A)	Noise Level(dBA)

- * Real-time display of parameters operation
- * With special designed rotors for 4*12 gel micro column blood cards.
- * Widely used in blood serology, blood type routine test, erythrocyte washing, micro column agglutination, immunology and other tests.
- * TFT Touch Screen

Intuitive control interface with numeric input buttons allows easy setting and editing of running conditions.

Technical Specifications:

Control	Microprocessor control	Program	30
Motor	Inverter motor direct drive	ACC/DEC	1-10 rates
Power supply	AC220V,50HZ,15A,1050W	Weight without rotor	135kg
Dimension(mm)	720*620*890mm(LxWxH)		

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Swing rotor	4×12 cards	4000r/min	2773xg	12 cards / box for swing bucket