

Benchtop • Low Speed • Ventilated • Centrifuge

T4-S



Scan to watch product video



4000 rpm	Maximum Speed
2100 xg	Maximum RCF
8x15 ml	Maximum Capacity
≤55dB (A)	Noise Level(dBA)

Technical Specifications:

Time range	0-99minutes adjustable	Speed accuracy	±10rpm
Power supply	AC220V,50HZ,5A,60W	Weight without rotor (kg)	15kg
Dimension	310x270x205mm(LxWxH)		

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Angle rotor	8×15/10/5ml	4000r/min	1900xg	Blood collection tube/PP/PC tube
	Adapter	8×2ml			Blood collection tube
NO.2 	Adapter	12×5ml	4000r/min	2100xg	Blood collection tube
	Adapter	12×2ml			Blood collection tube