

# PRP&PRF&CGF · Centrifuge

TPRP4-M



Scan to watch product video




|           |                  |
|-----------|------------------|
| 6000 rpm  | Maximum Speed    |
| 6080 xg   | Maximum RCF      |
| 8x100 ml  | Maximum Capacity |
| ≤58dB (A) | Noise Level(dBA) |

- \* Store up to 25 programs, easily switch between RPR, CGF, APRF, and IPRF protocols.
- \* Versatile variable-speed setting, compatible with 60ml/20ml/10ml/8ml injection syringes or PRP separation.

## Technical Specifications:

|               |                            |                            |            |
|---------------|----------------------------|----------------------------|------------|
| Time range    | 1-99minutes 59s adjustable | Speed accuracy             | ±10rpm     |
| Power supply  | AC220V,50HZ,10A, 570W      | ACC/DEC                    | 1-10 rates |
| Dimension(mm) | 580x440x360mm(LxWxH)       | Net weight (without rotor) | 40kg       |





Matched Rotors

| No.                                                                                                            | Rotor       | Capacity          | Max.speed | Max. RCF | Note                          |
|----------------------------------------------------------------------------------------------------------------|-------------|-------------------|-----------|----------|-------------------------------|
| <div>NO.1</div> <div></div>   | Angle rotor | 12x15/10/5ml      | 6000r/min | 6080xg   | PP/PC tube                    |
| <div>NO.2</div> <div></div> | Swing rotor | steel frame rotor | 4000r/min | 3041xg   | PP/PC tube                    |
|                                                                                                                | Swing frame | 4x100ml           |           |          | PP/PC tube                    |
|                                                                                                                | Adapter     | 4x50ml            |           |          | PP/PC tube                    |
|                                                                                                                | Adapter     | 4x30/20ml         |           |          | PP/PC tube                    |
|                                                                                                                | Adapter     | 4x8/10ml          |           |          | PP/PC tube                    |
| <div>NO.3</div> <div></div> | Swing rotor | steel frame rotor | 4000r/min | 3041xg   | Injection syringe or PRP tube |
|                                                                                                                | Swing frame | 4x2x100ml         |           |          | PP/PC tube                    |
|                                                                                                                | Adapter     | 4x2x50ml          |           |          | PP/PC tube                    |
|                                                                                                                | Adapter     | 4x2x30/20ml       |           |          | PP/PC tube                    |
|                                                                                                                | Adapter     | 4x2x8/10ml        |           |          | PP/PC tube                    |
| <div>NO.4</div> <div></div> | Angle Rotor | 8x8/10ml          | 3300r/min | 1339xg   | CGF special rotor and program |