

Benchtop • Low Speed • Refrigerated • Centrifuge

TR6-L



Scan to watch product video



6000 rpm	Maximum Speed
6080 xg	Maximum RCF
4x750 ml	Maximum Capacity
≤58dB (A)	Noise Level(dBA)

- * Temperature up to -20℃ Celslus
- * Accuracy of ±1℃
- * Fixed angle rotor, Swing out rotor, Microplate rotor

Technical Specifications:

Time range	1-99minutes 59s adjustable	Speed accuracy	±10rpm
Power supply	AC220V,50HZ,10A,1300W	Acceleration/Deceleration rates	1-10
Dimension	670x740x370mm(LxWxH)	Net weight(without rotor)	100kg

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
 NO.1	Angle rotor	12x15/10/5ml	6000r/min	6080xg	PP/PC tube
	Adapter	12x2ml			PP/PC tube
 NO.2	Swing rotor	steel frame rotor	5000r/min	4752xg	PP/PC tube
	Swing frame	4x100ml			PP/PC tube
	Swing frame	4x50ml			PP/PC tube
	Adapter	4x15/10ml			PP/PC tube

NO.3 	Swing rotor	steel frame rotor	4000r/min	3041xg	PP/PC tube
	Swing frame	4x2x100ml			PP/PC tube
	Swing frame	4x2x50ml			PP/PC tube
	Adapter	4x2x15/10ml			PP/PC tube
	Swing frame	4x8x15/10/5ml			PP/PC tube
	Rack	4×12×5/2ml			PP/PC tube
NO.4 	Swing bucket	72x5/2ml	4000r/min	2773xg	Vacuum blood collection tubes
NO.5 	Swing bucket	4×600ml	4000r/min	3354xg	Square cup
	Adapter	4×25×5/2ml			Vacuum blood collection tubes
	Adapter	4×12×15/10ml			PP/PC tube
	Adapter	4×5×50ml			PP/PC tube
	Adapter	4×4×100ml			PP/PC tube
NO.6 	Swing bucket	4×500ml	4000r/min	3354xg	Round cup
	Adapter	4×17×5/2ml			PP/PC tube
	Adapter	4×8×15/10ml			PP/PC tube
	Adapter	4×3×50ml			PP/PC tube
	Adapter	4×100ml			PP/PC tube
NO.7 	Swing bucket	4×750ml	4000r/min	3443xg	Round cup
	Adapter	4×26×5/2ml			PP/PC tube
	Adapter	4×13×15/10ml			PP/PC tube
	Adapter	4×5×50ml			PP/PC tube
	Adapter	4×3×100ml			PP/PC tube
NO.8 	Microplate rotor	2×2×96hole	4000r/min	2859xg	Square cup
	Adapter	2×32×1.5ml			PP/PC tube