

Benchtop • Low Speed • Ventilated • Centrifuge

T4-M



Scan to watch product video



4000 rpm	Maximum Speed
2680 xg	Maximum RCF
Angle rotor: 12x20 ml	Maximum Capacity
≤55dB (A)	Noise Level(dBA)

Technical Specifications:

Time range	0-99minutes adjustable	Speed accuracy	±10rpm
Power supply	AC220V,50HZ,5A,60W	Weight without rotor	18kg
Dimension	430x320x250mm(LxWxH)		

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Angle rotor	12×10/5ml	4000r/min	2680xg	PP/PC tube/ Vacuum blood collection tubes
	Adapter	12×2ml			Vacuum blood collection tubes
NO.2 	Angle rotor	12×20/15/10/5ml	4000r/min	2680xg	Vacuum blood collection tubes
	Adapter	12×2ml			Vacuum blood collection tubes
NO.3 	Angle rotor	18×10/5ml	4000r/min	2680xg	PP/PC tube/ Vacuum blood collection tubes
	Adapter	18×2ml			Vacuum blood collection tubes
	Adapter	18×1.5ml			PP/PC tube
NO.4 	Angle rotor	4×50ml	4000r/min	2680xg	PP/PC tube
	Adapter	4×15/10ml			Vacuum blood collection tubes
	Adapter	4×5/2ml			Vacuum blood collection tubes