

Benchtop • High Speed • Centrifuge

T21-MU



Scan to watch product video




21000 rpm	Maximum Speed
31554 xg	Maximum RCF
4x750ml	Maximum Capacity
≤58dB(A)	Noise(dBA)

Technical Specifications:

Time range	0-99minutes 59s adjustable	Speed accuracy	±10r/min
Motor	Motor Frequency Converter motor	Acceleration/Deceleration rates	1-10
Power supply	AC220V,50HZ,10A,770W	Display	LCD
Dimension	560x440x360mm(LxWxH)	Net weight(without rotor)	50kg

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
 NO.1	Angle rotor	12×1.5/2.2ml	21000r/min	31061xg	PP/PC tube
	Adapter	12×0.5/0.2ml			PP/PC tube
 NO.2	Angle rotor	12×5ml	21000r/min	31554xg	PP/PC tube
	Adapter	12×1.5ml			PP/PC tube
 NO.3	Angle rotor	24×1.5/2.2ml	16000r/min	23800xg	PP/PC tube
	Adapter	24×0.5/0.2ml			PP/PC tube
 NO.4	Angle rotor	10×15ml	16000r/min	31197xg	PP/PC tube
	Adapter	10×10ml			PP/PC tube
	Adapter	10×5ml			PP/PC tube
	Adapter	10×1.5ml			PP/PC tube
 NO.5	Angle rotor	48×1.5/2.2ml	15000r/min	23394xg	PP/PC tube
	Adapter	48×0.5/0.2ml			PP/PC tube
 NO.6	Angle rotor	12×10ml	15000r/min	20557xg	PP/PC tube
	Adapter	12×5ml			PP/PC tube

NO.7 	Angle rotor	6×50ml	15000r/min	22640xg	PP/PC tube
	Adapter	6×15ml			PP/PC tube
	Adapter	6×10ml			PP/PC tube
	Adapter	6×1.5ml			PP/PC tube
NO.8 	Angle rotor	6×100ml	13000r/min	20406xg	PP/PC tube
	Adapter	6×50ml			PP/PC tube
	Adapter	6×15ml			PP/PC tube
	Adapter	6×10ml			PP/PC tube
	Adapter	6×5/1.5ml			PP/PC tube
NO.9 	Angle rotor	8×50ml	13000r/min	20122xg	PP/PC tube
	Adapter	8×15ml			PP/PC tube
	Adapter	8×10ml			PP/PC tube
	Adapter	8×1.5ml			PP/PC tube
NO.10 	Swing bucket	4×250ml	4000r/min	3049xg	Round cup
	Adapter	4×9×5/2ml			Vacuum blood collection tube
	Adapter	4×5×15/10ml			PP/PC tube
	Adapter	4×50ml			PP/PC tube
	Adapter	4×100ml			PP/PC tube
	Adapter	4×9×1.5ml			PP/PC tube
NO.11 	Swing bucket	4×400ml	4000r/min	3049xg	Square cup
	Adapter	4×16×5/2ml			Vacuum blood collection tube
	Adapter	4×9×15/10ml			PP/PC tube
	Adapter	4×4×50ml			PP/PC tube
	Adapter	4×100ml			PP/PC tube
	Adapter	4×15×1.5ml			PP/PC tube
NO.12 	Swing bucket	4×500ml	4000r/min	3354xg	Round cup
	Adapter	4×17×5/2ml			Vacuum blood collection tube
	Adapter	4×8×15/10ml			PP/PC tube
	Adapter	4×3×50ml			PP/PC tube
	Adapter	4×100ml			PP/PC tube
NO.13 	Swing bucket	4×600ml	4000r/min	3354xg	Square cup
	Adapter	4×25×5/2ml			Vacuum blood collection tube
	Adapter	4×12×15/10ml			PP/PC tube
	Adapter	4×5×50ml			PP/PC tube
	Adapter	4×4×100ml			PP/PC tube
NO.14 	Swing bucket	4×750ml	4000r/min	3443xg	Round cup
	Adapter	4×26×5/2ml			Vacuum blood collection tube
	Adapter	4×13×15/10ml			PP/PC tube
	Adapter	4×5×50ml			PP/PC tube
	Adapter	4×3×100ml			PP/PC tube
NO.15 	Microplate rotor	2×2×96hole	4200r/min	3057xg	Square cup
	Adapter	2×32×1.5ml			PP/PC tube