


Benchtop • High Speed • Refrigerated • Centrifuge

TR21HI-M



Scan to watch product video



21000 rpm	Maximum Speed
31061 xg	Maximum RCF
4x100 ml	Maximum Capacity
≤58dB (A)	Noise Level(dBA)
-20℃~40℃	Temperature up to
±1℃	Celsius, Accuracy of


- * TFT Touch Screen: Intuitive control interface with numeric input buttons allows easy setting and editing of running conditions.
- * RTC (Real-Time Control) Function

Technical Specifications:

Time range	1-99minutes 59s adjustable	Speed accuracy	±10r/min
Power supply	AC220V,50HZ,10A, 650W	ACC/DEC	1-10 rates
Dimension	560x340x330mm(LxWxH)	Net weight(without rotor)	42kg

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
	Angle rotor	12×1.5/2.2ml	21000r/min	31061xg	PP/PC tube
	Adapter	12×0.5ml			PP/PC tube
	Adapter	12×0.2ml			PP/PC tube

NO.2 	Angle rotor	12x5ml	16500r/min	19480xg	PP/PC tube
	Adapter	12x1.5ml			PP/PC tube
NO.3 	Angle rotor	24x1.5/2.2ml	15000r/min	20124xg	PP/PC tube
	Adapter	24x0.5ml			PP/PC tube
	Adapter	24x0.2ml			PP/PC tube
NO.4 	Angle rotor	12x8x0.2ml	14000r/min	23120xg	PCR8 joint pipe
NO.5 	Angle rotor	12×10ml	15000r/min	20557xg	PP/PC tube
	Adapter	12x5ml			PP/PC tube
	Adapter	12x1.5ml			PP/PC tube
NO.6 	Angle rotor	6×50ml	12000r/min	14489xg	PP/PC tube
	Adapter	6×15ml			PP/PC tube
	Adapter	6x10ml			PP/PC tube
	Adapter	6x1.5ml			PP/PC tube
NO.7 	Angle rotor	8×15ml	11000r/min	13528xg	PP/PC tube
	Adapter	8×5ml			PP/PC tube
	Adapter	8×1.5ml			PP/PC tube
NO.8 	Angle rotor	4×10/100ml	11000r/min	13280xg	PP/PC tube
	Adapter	4×50ml			PP/PC tube
	Adapter	4×15ml			PP/PC tube
	Adapter	4×10ml			PP/PC tube
	Adapter	4×5ml			PP/PC tube
	Adapter	4×1.5ml			PP/PC tube