

Gel Card Centrifuge

XK-12B



Scan to watch product video



3500 rpm	Maximum Speed
1610xg	Maximum RCF
12 gel cards	Capacity
≤50dB (A)	Noise Level(dBA)

- * Real-time display of multiple parameters operation.
- * Automatic door opening alarm after shutdown.
- * Professional programmed design, no manual setup required, simply run it directly.
- * Widely used in blood serology, blood type routine test, erythrocyte washing, micro column agglutination, immunology and other tests.

Technical Specifications:

Time Range	0-99min	Max Capacity	12 cards
Motor	Frequency conversion A.C. motor	Program	±10r/min
Power supply	AC220V,50HZ,5A,60W	Display	LCD
Dimension	430*320*250mm	Weight without rotor	13kg

Matched Rotors

No.	Item	Capacity	Max.speed	Max. RCF	Time	Note
NO.1 	Pre-centrifuge program group	12 cards	3500r/min	1232xg	1 minute	Program group can be set for different blood type cards.
NO.2 	Test program group	12 cards	900r/min	81xg	2 minutes	
			1500r/min	226xg	3 minutes	