

Benchtop • High Speed • Ventilated • Centrifuge

MI16



Scan to watch product video



16500 rpm	Maximum Speed
19480 xg	Maximum RCF
12x5 ml	Maximum Capacity
≤55dB (A)	Noise Level(dBA)

Technical Specifications:

Time range	1-99minutes 59s adjustable	Speed accuracy	±10r/min
Power supply	AC220V,50HZ,10A, 370W	ACC/DEC	1-10 rates
Dimension	340x260x220mm(LxWxH)	Weight without rotor	15kg
Max. acceleration time [s]	15 seconds	Max. braking time [s]	15 seconds

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Angle rotor	12×1.5/2.2ml	16500r/min	19175xg	EP/PP tube
	Adapter	12×0.5/0.2ml			EP/PP tube
NO.2 	Angle rotor	12×5ml	16500r/min	19480xg	EP/PP tube
	Adapter	12×1.5ml			EP/PP tube
NO.3 	Angle rotor	8×8×0.2ml	13500r/min	16480xg	PCR 8 strip tubes
NO.4 	Angle rotor	24×1.5-2.2ml	14000r/min	17530xg	EP/PP tube
	Adapter	24×0.5/0.2ml			EP/PP tube