

Floor Standing • Low Speed • Refrigerated • Centriuge

FR8-B12



Scan to watch product video



8000 rpm	Maximum Speed
14900 xg	Maximum RCF
6x2400 ml	Maximum Capacity
12 pcs 500ml blood bags	
≤65dB (A)	Noise Level(dBA)
-20℃~40℃	Temp.range Accuracy of ±1℃

- * RTC (Real-Time Control) Function
- * Castor wheels for ease of mobility

Technical Specifications:

Time range	1-99h59minutes adjustable	Acceleration/Deceleration rates	1-- 10
Dimension	940×890×1000mm(L×W×H)	Daily-use program	30
Display	LCD	Net weight(without rotor)	600kg
Speed accuracy	±20r/min	Power supply	AC380V,50HZ,15A
Control&drive Microprocessor	Brushless A.C. motor with big torque and frequency conversion,microprocessor control		

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Angle Rotor	6×1000ml	8000r/min	14900xg	
NO.2 	Swing Rotor	6×2400ml	4000r/min	5330xg	12 pcs of 500ml blood bags
NO.3 	Swing Rotor	6×2×1000ml	4000r/min	5800xg	