

Benchtop • Low Speed • Ventilated • Centrifuge

T5-L



Scan to watch product video



5000 rpm	Maximum Speed
3354xg	Maximum RCF
4x100 ml	Maximum Capacity
≤55dB (A)	Noise Level(dBA)

Technical Specifications:

Time range	0-99minutes 59s	Speed accuracy	±10 rpm
Motor	Frequency Converter motor	Acceleration / Deceleration rates	0-9
Power supply	AC220V,50HZ,5A,320W	Display	LCD
Dimension	500x400x300mm(LxWxH)	Weight without rotor	25kg



Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Angle rotor	12×15/10/5ml	5000r/min	3047xg	PP/PC tube
	Adapter	12×2ml			Vacuum blood collection tube
NO.2 	Angle rotor	6×50ml	5000r/min	3354xg	PP/PC tube
	Adapter	6×15/10ml			PP/PC tube
	Adapter	6×5ml			Vacuum blood collection tube
NO.5 	Swing rotor	Steel frame rotor	4000r/min	2593xg	PP/PC tube
	Frame	4×100ml			PP/PC tube
	Frame	4×50ml			PP/PC tube
	Adapter	4×10-15ml			Vacuum blood collection tube
	Adapter	4×4×5-10-15ml			Vacuum blood collection tube
	Adapter	4×6×5-10ml			PP/PC tube