

Bio-safety · Centrifuge

TBSA6-L



Scan to watch product video





6000 rpm	Maximum Speed
6080 xg	Maximum RCF
4x600 ml	Maximum Capacity
≤58dB (A)	Noise Level(dBA)

- * TFT Touch Screen
- Intuitive control interface with numeric input buttons allows easy setting and editing of running conditions.
- * Gas-tight head for high-temp sterilization
- * UV function in the centrifugal cavity
- * Instant disinfection at startup

Technical Specifications:

Time range	0-99minutes 59s adjustable	ACC/DEC	0-9
Speed accuracy	±10r/min	Net weight(without rotor)	50kg
Power supply	AC220V,50HZ,10A, 770W	Dimension	600x470x370mm(LxWxH)



Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Angle rotor	12x15/10/5ml	6000r/min	6080xg	PP/PC tube
	Adapter	12x2ml			
NO.2 	Swing bucket	4x600ml	4000r/min	3354xg	Square bucket
	Adapter	4x25x5/2ml			
	Adapter	4x12x15/10ml			
	Adapter	4x5x50ml			
	Adapter	4x4x100ml			
NO.3 	Swing bucket	4×500ml	4000r/min	3354xg	Round bucket with air tight cover
	Adapter	4×24×5/2ml			
	Adapter	4×13×15/10ml			
	Adapter	4×4×50ml			
	Adapter	4×2×100ml			