

Benchtop • High Speed • Ventilated • Centrifuge

MI16H



Scan to watch product video




16500 rpm	Maximum Speed
19480 xg	Maximum RCF
10x5 ml	Maximum Capacity
≤55dB (A)	Noise Level(dBA)

Technical Specifications:

Time range	1-99minutes 59s adjustable	Speed accuracy	±10r/min
Power supply	AC220V,50HZ,10A, 400W	ACC/DEC	1-10 rates
Dimension	340*260*220mm(LxWxH)	Weight without rotor	15kg

Matched Rotors

No.	Rotor	Capacity	Max.speed	Max. RCF	Note
NO.1 	Angle rotor	12×1.5/2.2ml	16500r/min	19175xg	PP/PC tube
	Adapter	12×0.5/0.2ml			PP/PC tube
NO.2 	Angle rotor	10×5ml	16500r/min	19480xg	PP/PC tube
	Adapter	10×1.5ml			PP/PC tube
NO.3 	Capillary rotor	24 pieces	12000r/min	13684xg	Capillary tube